

European Spine Phantom

QRM-ESP

The standard for quality control in bone mineral densitometry by DXA and QCT.

The European Spine Phantom QRM-ESP is now the standard for quality control - accuracy and reproducibility - in spinal BMD quantification for more than a decade [3]. Its design is closely anthropomorphic by its outer measures and internal lumbar vertebrae inserts such that standard patient protocols can be utilized, both in DXA and QCT.

The phantoms main body consists of waterequivalent resin, the three inserts contain varying amounts of calcium hydroxyapatite (CaHA) to cover the full physiological range of spongeous and cortical bone densities for all age groups.

QRM-ESP allows to check reproducibility and accuracy with the following quantities:


- ✓ bone mineral content (BMC) in g for DXA
- √ bone mineral area density (BMD) in g/cm2
- ✓ for DXA
- ✓ projected area (A) in cm2 for DXA
- trabecular and cortical bone mineral density in g/cm3 for QCT
- ✓ cortical thickness in mm for QCT
- ✓ positioning accuracy in QCT

Cross-sectional view of the three vertebrae


L1 Spongious: 50 HA (mg/cm³)
Cortical: 800 HA (mg/cm³)
L2 Spongious: 100 HA (mg/cm³)
Cortical: 800 HA (mg/cm³)
L3 Spongious: 200 HA (mg/cm³)
Cortical: 800 HA (mg/cm³)



Resin based water-equivalent phantom body approximately oval shape but with flattend sides for ease of positioning.




DXA scan of QRM-ESP



Description see left side



European Spine Phantom



CT images of the QRM-ESP. Scans of the midvertebral selections with automatically determined ROI for spongious and cortical bone superimposed.



DXA scan of QRM-ESP

Specification

Area density of vertebrae (AP):

0.5, 1.0, 1.5 (g/cm²)

HA in cortical structures and spinal process:

400, 800 HA (mg/cm³)

Accuracy $\pm 3\%$ of specified values $\pm 1\%$ of certified values

References

- [1] Kalender, W.A.: A Phantom for Standardization and Quality Control in Spinal Bone Mineral Measurements by QCT and DXA. Design Considerations and Specifications. Med. Phys. 19 (1991) 583-588
- [2] Genant, H.K. et al.: Universal Standardization for Dual X-Ray Absorptiometry: Patient and Phantom Cross-Calibration Results. J.Bone Miner. Res. 9 (1994) 1503-1514
- [3] Kalender, W.A., Felsenberg, D., Genant, H., Fischer, M., Dequeker, J., Reeve; J.: The European Spine Phantom a tool for standardization and quality control in spinal bone measurements by DXA and QCT. European J. Radiology 20 (1995) 83-92